***Storage System***  
**Ultrasound Probe Cabinet**

**Operator's Manual**

**Manufactured by:**

**CIVCO Medical Solutions**  
750 Lois Lane  
Coralville, IA 52241 USA  
319.248.6757 / 800.445.6741  
[WWW.CIVCO.COM](http://WWW.CIVCO.COM)

Copyright © 2021 All Rights Reserved.  
GUS® is a registered trademark of CIVCO Medical Solutions.  
All other trademarks are property of their respective owners.  
Printed in the USA.

# Table of Contents

<b>General Information</b> .....	<b>2</b>
<i>Intended Use</i> .....	2
<i>User Qualifications</i> .....	3
<i>Warning Symbols</i> .....	3
<i>Signal Words</i> .....	3
<b>Installation and Setup</b> .....	<b>4</b>
<i>Upon Delivery</i> .....	4
<i>Wall Mounting Instructions</i> .....	4
<b>Instructions for Use</b> .....	<b>5</b>
<i>Operation</i> .....	5
<i>Probe Hanging: Vaginal, Rectal &amp; General Purpose Probes Only</i> .....	6
<i>Probe Hanging: TEE Probes Only</i> .....	7
<b>Cleaning the System</b> .....	<b>7</b>
<b>Routine Maintenance</b> .....	<b>8</b>
<i>Daily System Inspection</i> .....	8
<i>Filter Maintenance</i> .....	8
<b>Limited Warranty</b> .....	<b>9</b>
<b>Technical Support</b> .....	<b>9</b>
<b>Registration Card</b> .....	<b>9</b>
<b>Appendix A: Troubleshooting</b> .....	<b>10</b>
<b>Appendix B: Site Requirements</b> .....	<b>11</b>
<i>Actual Size</i> .....	11
<i>Minimum Space Required</i> .....	11
<i>Electrical Requirements</i> .....	11
<i>Environmental Conditions</i> .....	11
<b>Appendix C: Product Specifications</b> .....	<b>12</b>
<i>Product Consumables</i> .....	12

## ***General Information***

This operator's manual contains essential information regarding the proper use of this system. Before use, thoroughly review this manual. Use all equipment as instructed.

Keep this manual and all related instruction manuals in a safe and accessible location.

If you have any questions or comments regarding any information in this operator's manual, please contact CIVCO Medical Solutions.

## ***Intended Use***

The CIVCO Storage System provides for a safe and clean storage environment for your disinfected ultrasound probes (See Figures 1 and 2). A blower draws air through a replaceable HEPA filter providing a continuous flow of clean filtered air allowing safe and clean storage. The easily accessible filter is located on top of the Storage System. The Storage System's locked door provides secure and safe storage. The filter life is monitored and the indicator light will signal when a replacement filter is required.

**WARNING: Do not use this system for any purpose other than its intended use.**

Figure 1: Vaginal, Rectal and General Purpose Probe Storage System

Figure 2: TEE Probe Storage System

## ***User Qualifications***

The CIVCO Storage System is intended for use in non-sterile healthcare settings by trained, qualified healthcare personnel.

## ***Warning Symbols***

The meaning of the labels and symbols that appear on the packaging and/or the system are as follows:

Dangerous electrical voltage

Hazardous area

Protective Grounding

## ***Signal Words***

The following signal words are used throughout this manual:

**WARNING:** Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

**CAUTION:** Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices or potential equipment damage.

**NOTE:** Indicates additional helpful information.

### **WARNINGS:**

**WARNING:** Do not use this system until it has been properly installed.

**WARNING:** Do not operate this device if it is not working properly or if any part of the device has been dropped or damaged.

**WARNING:** Connect the system to a hospital grade GFCI receptacle only.

**WARNING:** Follow the probe manufacturer's instructions for storage and handling.

### **CAUTIONS:**

**CAUTION:** Do not place any objects on top of the system. This will affect the airflow of the system.

---

## ***Installation and Setup***

### ***Upon Delivery***

Please check to ensure your new system arrived with the following parts:

#### **COMPONENT LIST:**

- Either:
  - Storage System for TEE probes with shelving for holding up to three (3) transesophageal ultrasound probes.
- Or:
  - Storage System for VR and general purpose probes with shelving for holding up to six (6) vaginal, rectal or general purpose ultrasound probes.
- Wall bracket with supporting wall anchors

### ***Wall Mounting Instructions***

**WARNING: Do not use this system until it has been properly installed.**

**WARNING: Do not operate this device if it is not working properly or if any part of the device has been dropped or damaged.**

**WARNING: Connect the system to a hospital grade GFCI receptacle only.**

**WARNING: To prevent electrical or wall damage, it is recommended the Storage System be installed by the facilities department.**

**CAUTION: Use two or more people to lift and position the system.**

1. You will need the following items before attempting to install this unit:

- Measuring tape
- Level
- Pencil or pen

**NOTE: The mounting assembly consists of a bracket pre-mounted to the CIVCO Storage System and a separate wall bracket to be mounted securely to the wall. The unit must be attached to the wall using the interlocking mounting bracket supplied prior to use.**

2. Locate a wall with adequate clearance for the cabinet. Ensure the wall has sufficient strength to mount the CIVCO Storage System. See "Appendix B: Site Requirements".
3. Make sure that there is a minimum of 6" clearance on the right side of the cabinet for the power cord and On/Off switch.
4. Mark the center point of the cabinet on the wall.

5. From the floor measure up 64" and make a mark on the floor.
6. From this point make a level horizontal line on the wall extending 11 ½" from this center point in both directions. The length of this horizontal line should be 23" total.
7. Hold the mounting bracket with its bottom edge on the horizontal line from step 6 and the center mounting hole aligned with the center mark of the cabinet. Mark the position of the mounting holes on the bracket.
8. Attach the bracket to the wall using the mounting hardware supplied at the points marked from step 7. Use EZ anchors (supplied) if no wall studs are available for the mounting.
9. Once the bracket has been securely mounted to the wall, lift up the device using the lift handles on the side of the device.
10. Push the device against the wall and let it down easily to engage the wall bracket.
11. Double check that the bracket is engaged.
12. Once the wall bracket is properly engaged, attach the power cord to the receptacle on the right side of the device.

**NOTE: The inlet power module which contains the main power switch, fuses and line cord receptacle are located on the upper right side of the device. The light in the power switch will be green when the system is on. The ratings plate and serial number are located on the right side of the system below the power switch receptacle.**

## ***Instructions for Use***

### ***Operation***

**WARNING: Disinfect system prior to use. See "*Cleaning the System*" section for instructions for cleaning.**

**WARNING: Keep system door closed while in use.**

**NOTE: Keep system on while probe(s) are being stored.**


1. Ensure the power cord is attached to the receptacle on the right side of the device.
2. Plug the power cord into a Hospital Grade receptacle and switch on the power switch located in the power control module above the power cord. The light in the power switch will be green when the power is on.
3. The blower will start when the power switch is turned on.
4. Unlock the system door with the key supplied and push the key button/latch to open the door. Make sure the key is kept in a safe place.
5. On the upper right of the system is a rectangular indicator light to notify operator when it is necessary to change the HEPA filter.
6. When filter is within operating parameters (11-months) the indicator light will be green.
7. Once eleven (11) months running time has elapsed, the indicator light will flash red indicating it is time to change the filter. See "*Filter Maintenance*" section for instructions on replacing the filter.

## ***Probe Hanging: Vaginal, Rectal & General Purpose Probes Only***

**WARNING:** Follow the probe manufacturer's instructions for storage and handling.

**CAUTION:** When installing the Vaginal/Rectal/General Purpose probes, ensure at all times that the connector housing and cable are kept away from the transducer to avoid cross-contamination.

1. Open the system door by depressing the key button/latch.
2. While holding the transducer in an upright position, guide it up through the circular slot in the upper shelf. Rotate the transducer distal tip/lens so that it is facing outward.
3. Guide the cable through the slots in the middle and lower shelves. Rest the transducer securely on the middle shelf.
4. Insert the connector housing into the circular slot on the lower shelf, with the cable looped and hanging in a downward position behind the connector.
5. Repeat steps 2 through 4, as needed.
6. Close door.

**NOTE:** Probe manufacturers will have the connector housing orientated in different positions. Ensure that the connector housing is secured safely on the lower shelf before moving to the next step.

### **FOR GENERAL PURPOSE PROBES:**

1. Open the system door by depressing the key button/latch.
2. While holding the transducer, guide the cable through the slots in the middle and lower shelves. Rest the transducer securely on the middle shelf.
3. Insert the connector housing into the circular slot on the lower shelf, with the cable looped and hanging in a downward position behind the connector.

**NOTE:** Probe manufacturers will have the connector housing orientated in different positions. Ensure that the connector housing is secured safely on the lower shelf.

4. Repeat Steps 2 and 3, as needed.


## Probe Hanging: TEE Probes Only

**WARNING:** Follow the probe manufacturer's instructions for storage and handling.

**CAUTION:** When installing the TEE probe, ensure at all times that the connector housing and cable are kept away from the shaft and transducer lens to avoid cross-contamination.

**CAUTION:** Ensure that the distal tip and transducer lens are stored securely away from the system door when closed.

1. Open the system door by depressing the key button/latch.
2. Hold the TEE Probe above the upper shelf while inserting the shaft into the slot on the middle shelf, resting the distal tip/transducer lens into the lower slot on the lower shelf.
3. Lower the control housing strain relief into the circular slot located on the upper shelf.
4. Take the cable and loop it over the strain relief located above the upper shelf.
5. Insert the cable into the circular slots located on the upper and middle shelf.
6. Guide the cable through the circular slot on the middle shelf, and place the connector housing on the upper shelf on the left side of the TEE probe.
7. Repeat Steps 2 through 6, as needed.
8. Close the door.

## Cleaning the System

**WARNING:** Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

**CAUTION:** Disconnect power when cleaning the system.

**CAUTION:** Do not use solvents or isopropyl alcohol to clean the acrylic door. Doing so may damage the doors. Use a mild detergent solution on a soft cleaning cloth.

**CAUTION:** Do not wipe surfaces with hard or abrasive materials. Doing so will scratch the surfaces.

1. To clean metal surfaces, use a surface disinfectant authorized by your infection control department.
2. Make sure that the device is completely dry before use.

---

## ***Routine Maintenance***

This system requires minimal routine maintenance and inspection.

### ***Daily System Inspection***

1. Visually inspect electrical connections prior to operation.
2. Check area around the device to ensure safe operation.
3. Check blower operation by placing hand over the outlet on the top of the inside of the cabinet and feeling for air movement.
4. Check blower operation by monitoring for unusual or harsh noises.
5. Confirm maintenance indicator is green; this indicator is located on the upper right side of the system behind the power inlet module. If the indicator is red, the HEPA filter needs to be replaced. See "*Filter Maintenance*" section for instructions on replacing the filter.

### ***Filter Maintenance***

Filter life is monitored by a timer; replacement is signaled after the filter has been used for eleven (11) months. The timer is only active when the blower is on. Call your CIVCO Medical representative or local distributor to order replacement filters.

**CAUTION: Filter life is dependent upon frequency and conditions of use. As a result, the filter may need to be replaced prior to the 11 month guideline.**

**CAUTION: Disconnect power when replacing the filter.**

1. Turn the main power switch off and disconnect the power cord from the electrical outlet.
2. Due to the height of the device, a stool or ladder is recommended to access the filter.
3. There are two latches on the top panel of the device. To open the top and access the filter, turn the latches 90 degrees. Use the latches on the top panel to lift the top panel to a vertical position.
4. Remove the used filter by pulling it up.
5. Put the used filter into a plastic bag and dispose of it with other non-hazardous waste. **This does not need to be red tagged.**
6. Install the new HEPA filter with the gasket and arrow pointing down.

**NOTE: HEPA filters are very fragile. When handling the HEPA filter, use the outer metal frame to grip the filter when changing or installing. Do not damage the filter material. Handle the filter on the outer metal frame.**

7. Once the filter is in place, swing the top panel down, push down and secure by closing the left and right latches previously opened.
8. Pull up on the latches on the top of the device to make sure latches are secure.

9. Plug the device back into the wall outlet and turn on the power switch.

**NOTE: The filter change indicator light will still be flashing red at this point and will need to be reset.**

10. To reset the filter change indicator light, located on the upper side of the system behind the power inlet module, hold the button down for 5 seconds until the maintenance indicator light turns green. This will reset the filter change indicator light for another 11 months.

### ***Limited Warranty***

CIVCO Medical Solutions warrants the CIVCO Storage System to be free of defects in materials and workmanship under normal use, for a period of one year from the date of delivery. This warranty does not include filters or acrylic door, unless it can be determined that the failure is due to defects in material or workmanship.

If the system is not working correctly or is defective, contact CIVCO Medical Solutions and we will repair or replace the device at our option during the warranty period.

CIVCO Medical Solutions assumes no liability for consequential damages of any kind as a result of the use or misuse of the system by the purchaser, the purchaser's employees, or any others.

**THIS WARRANTY IS EXPRESSLY IN LIEU OF ALL OTHER WARRANTIES, GUARANTEES, OBLIGATIONS OR OTHER LIABILITIES EXPRESSED OR IMPLIED, AND THERE ARE NO WARRANTIES WHICH EXTEND BEYOND THE DESCRIPTION ON THE FACE HEREOF. WARRANTIES OF MERCHANTABILITY AND OF FITNESS FOR A PARTICULAR PURPOSE ARISING FROM A COURSE OF DEALING OR USAGE ARE SPECIFICALLY EXCLUDED.**

This warranty does not cover conditions and damages resulting from any of the following:

- Improper installation, delivery within the facility, or maintenance.
- Any repair, modification, alteration or adjustment not authorized by the manufacturer.
- Misuse, abuse, accidents or unreasonable use.
- Incorrect electric current, voltage or supply.
- Failure to follow proper procedures for use outlined in this manual.
- Consequential or incidental damages sustained by any person as a result of any breach of this warranty.
- Some states do not allow the exclusion or limitation of consequential or incidental damages, so the above exclusion may not apply.

### ***Technical Support***

CIVCO Medical Solutions, 750 Lois Lane, Coralville, IA 52241 USA  
Tel: 319.248.6757 / 800.445.6741  
WWW.CIVCO.COM

### ***Registration Card***

Please complete and return the attached Warranty Card to CIVCO Medical Solutions. Call 1-800-445-6741 to purchase additional filters.

---

## ***Appendix A: Troubleshooting***

If the CIVCO Storage System blower is not running, check the following.

1. Ensure that the device is plugged into a live wall outlet.
2. Make sure the On/Off switch is on. The switch will be illuminated when on.
3. In the event that the system is still not working, you will need to check the fuses, which are externally mounted in the power entry module with the On/Off switch and power cord receptacle.
4. The fuses are located in a rectangular receptacle between the switch and the power cord. Remove the fuse receptacle by squeezing the two tabs on the fuse receptacle and pulling outward. If fuses are blown, replace fuse.
5. If the system is still not working, call CIVCO Medical Solutions at 800.445.6741 or call your local distributor.

**NOTE: If the filter change indicator light is flashing red, this indicates it is time to change the HEPA filter. Refer to "Filter Maintenance" section for instructions on replacing the filter.**

---

## ***Appendix B: Site Requirements***

### ***Actual Size***

- Width: 24.0" 60.96 cm
- Depth: 10.0" 25.40 cm
- Height: 74.0" 187.96 cm
- Weight: 100 lbs 45.36 kg

### ***Minimum Space Required***

- Width: 30.0" 76.20 cm
- Depth: 36.0" 91.44 cm
- Height: 84.0" 213.36 cm

### ***Electrical Requirements***

- Service: 110-120 VAC 60Hz, 15A Hospital Grade GFCI receptacle
- Fuse: 1.0 Amp 250V 5mm x 20mm FSF Panel Mounted
- Current Draw: 0.5 Amp

### ***Environmental Conditions***

- Indoor use only
- Altitude up to 2,000 meters
- Temperature range: 5°C to 40°C
- Maximum relative humidity: 80% for temperatures up to 31°C, decreasing linearly to 50% relative humidity at 40°C
- Main supply voltage fluctuations not to exceed +/- 10% of the nominal voltage
- Insulation over voltage category II
- Pollution degree 2

---

## ***Appendix C: Product Specifications***

### ***Product Consumables***

<b>CIVCO Part #</b>	<b>Description</b>
610-1328	Replacement Filter

**USA and Canadian Standards Approved.**

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