Storage System **Ultrasound Probe Cabinet**

Operator's Manual

Manufactured by:

CIVCO Medical Solutions
750 Lois Lane
Coralville, IA 52241 USA
319.248.6757
WWW.CIVCO.COM

Copyright © 2025 All Rights Reserved.
All other trademarks are property of their respective owners.
Printed in the USA.

Table of Contents

General Information	2
Intended Use	2
User Qualifications	3
Symbol Glossary	3
Signal Words	4
Installation and Setup	5
Upon Delivery	5
Wall Mounting Instructions	5
Reversing Cabinet Door	6
Instructions for Use	7
Operation	7
Probe Hanging: Vaginal, Rectal & General Purpose Probes Only	8
Probe Hanging: TEE Probes Only	9
Cleaning and Disinfecting the System	9
Routine Maintenance	10
Daily System Inspection	10
Filter Maintenance	10
Technical Support	11
Registration Card	11
Appendix A: Troubleshooting	12
Appendix B: Site Requirements	13
Actual Size	13
Minimum Space Required	13
Electrical Requirements	13
Environmental Conditions	13
Mounting and Installation Diagrams	14
Appendix C: Product Specifications	15
Product Consumables	15
Airflow Specifications	15
Noise Level Specifications	15
Cabinet Materials	15
Shelving Specifications	16
Appendix D: EMC Compliance Information	18

General Information

This operator's manual contains essential information regarding the proper use of this system. Before use, thoroughly review this manual. Use all equipment as instructed.

Keep this manual and all related instruction manuals in a safe and accessible location.

If you have any questions or comments regarding any information in this operator's manual, please contact CIVCO Medical Solutions.

Intended Use

The CIVCO Storage System provides for a safe and clean storage environment for your disinfected ultrasound probes (See Figures 1 and 2). A fan draws air through a replaceable HEPA filter providing a continuous flow of clean filtered air allowing safe and clean storage. The easily accessible filter is located on top of the Storage System. The Storage System's locked door provides secure and safe storage. The filter life is monitored and the indicator light will signal when a replacement filter is required.

WARNING: Do not use this system for any purpose other than its intended use.

Figure 1: General Purpose & Endocavity Ultrasound Probe Cabinet

Figure 2: TEE Ultrasound Probe Cabinet

User Qualifications

The CIVCO Storage System is intended for use in non-sterile healthcare settings by trained, qualified healthcare personnel.

Symbol Glossary

The meaning of the labels and symbols that appear on the packaging and/or the system are as follows:

Symbol	Title of Symbol	Description of Symbol
	CAUTION (ISO 15223-1, 5.4.4)	Indicates the need for user to consult instructions for use for important cautionary information such as warning and precautions that cannot, for a variety of reasons, be presented on the device itself.
	Manufacturer (ISO 15223-1, 5.1.1)	Indicates the device manufacturer.
	Date of manufacture (ISO 15223-1, 5.1.3)	Indicates the date when device was manufactured.
	Consult instructions for use (ISO 15223-1, 5.4.3)	Indicates the need for the user to consult instructions for use.
	Catalog number (ISO 15223-1, 5.1.6)	Indicates the manufacturer's catalog number so the device can be identified.
	Quantity (IEC 60878, 2794)	Indicates quantity of contents.
	Serial number (ISO 15223-1, 5.1.7)	Indicates manufacturer's serial number so a specific device can be identified.
	Batch code (ISO 15223-1, 5.1.5)	Indicates the manufacturer's batch code so the batch or lot can be identified.
	Dangerous Voltage (IEC 60417-5036)	Indicates hazards arising from dangerous voltages.
	Functional Grounding (IEC 60417-5018)	Identifies a functional grounding terminal, for example, of a specially designed grounding system to avoid causing malfunction of the equipment.
	NRTL/MET-C Product Certification	Indicates the product meets the minimum requirements of applicable safety standards.

Signal Words

The following signal words are used throughout this manual:

WARNING: Indicates a potentially hazardous situation which, if not avoided, could result in death or serious injury.

CAUTION: Indicates a potentially hazardous situation which, if not avoided, may result in minor or moderate injury. It may also be used to alert against unsafe practices or potential equipment damage.

NOTE: Indicates additional helpful information.

WARNINGS:

WARNING: Do not use this system until it has been properly installed.

WARNING: Do not operate this device if it is not working properly or if any part of the device is dropped or damaged.

WARNING: If system loses power unexpectedly, clean and disinfect both the system and devices prior to reuse.

WARNING: Connect the system to a hospital grade GFCI receptacle only.

WARNING: Follow the probe manufacturer's instructions for storage and handling.

WARNING: Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 12 inches (30 cm) to any part of the system, including cables specified by the manufacturer. Otherwise, degradation of the performance of this equipment could result.

WARNING: Do not operate the device if the fan or filter change indicator does not function or display correctly due to EM disturbances, as this may affect device performance and safety.

CAUTIONS:

CAUTION: Do not place any objects on top of the system. This will affect the airflow of the system.

CAUTION: Avoid operation adjacent to or stacked with other electrical equipment.

CAUTION: Use of the supplied power cord is required to maintain proper grounding and safe operation of the device. Use of power cords not supplied by CIVCO Medical Solutions may affect the safety and performance of this device.

NOTE: The emissions characteristics of this equipment make it suitable for use in industrial areas and hospitals (CISPR 11 class A). If it is used in a residential environment (for which CISPR 11 class B is normally required), this equipment might not offer adequate protection to radio-frequency communication services. The user might need to take mitigation measure, such as relocating or re-orienting the equipment.

Installation and Setup

Upon Delivery

Please check to ensure your new system arrived with the following parts:

COMPONENT LIST:

- Either:
 - Storage System for TEE probes with shelving for holding up to three (3) transesophageal ultrasound probes.
- Or:
 - Storage System for VR and general purpose probes with shelving for holding up to six (6) vaginal, rectal or general purpose ultrasound probes.
- Wall bracket with supporting wall anchors

Wall Mounting Instructions

FOR MOUNTING AND INSTALLATION DIAGRAMS, REFER TO APPENDIX B.

WARNING: Do not use this system until it has been properly installed.

WARNING: Do not operate this device if it is not working properly or if any part of the device is dropped or damaged.

WARNING: Connect the system to a hospital grade GFCI receptacle only.

WARNING: To prevent electrical or wall damage, it is recommended the Storage System be installed by the facilities department.

CAUTION: Use two or more people to lift and position the system.

1. You will need the following items before attempting to install this unit:

- Measuring tape
- Level
- Pencil or pen

NOTE: The mounting assembly consists of a bracket pre-mounted to the CIVCO Storage System and a separate wall bracket to be mounted securely to the wall. The unit must be attached to the wall using the interlocking mounting bracket supplied prior to use.

2. Locate a wall with adequate clearance for the cabinet. Ensure the wall has sufficient strength to mount the CIVCO Storage System. See "Appendix B: Site Requirements".
3. Make sure that there is a minimum of 6 inches (15.2 cm) clearance on the right side of the cabinet for the power cord and On/Off switch.

4. Mark the center point of the cabinet on the wall.
5. From the floor measure up 64 inches (162.6 cm) and make a mark on the floor.
6. From this point make a level horizontal line on the wall extending 11.5 inches (29.2 cm) from this center point in both directions. The length of this horizontal line should be 23 inches (58.4 cm) total.
7. Hold the mounting bracket with its bottom edge on the horizontal line from step 6 and the center mounting hole aligned with the center mark of the cabinet. Mark the position of the mounting holes on the bracket.
8. Attach the bracket to the wall using the mounting hardware supplied at the points marked from step 7. Use EZ anchors (supplied) if no wall studs are available for the mounting.
9. Once the bracket has been securely mounted to the wall, lift up the device using the lift handles on the side of the device.
10. Push the device against the wall and let it down easily to engage the wall bracket.
11. Double check that the bracket is engaged.
12. Once the wall bracket is properly engaged, attach the power cord to the receptacle on the right side of the device.

NOTE: The inlet power module which contains the main power switch, fuses and line cord receptacle are located on the upper right side of the device. The light in the power switch will be green when the system is on. The ratings plate and serial number are located on the right side of the system below the power switch receptacle.

Reversing Cabinet Door

1. Ensure cabinet is securely mounted to the wall.
2. Open system door by depressing key button/latch.
3. Locate spring latch pin on inside bottom of door and pull to compress.
4. With latch pin compressed, tilt door and pull to release from cabinet.
5. Invert door so spring latch pin faces up.
6. Align pin on bottom of door with corresponding hole on cabinet and insert.
7. Compress spring latch on top of door, align pin and release to secure door to cabinet.
8. Close door and ensure proper installation.

Instructions for Use

Operation

WARNING: Clean and disinfect system prior to use. See "*Cleaning and Disinfecting the System*" section for instructions for cleaning and disinfecting the system.

WARNING: Keep system door closed while in use.

NOTE: Keep system on while probe(s) are being stored.

1. Ensure the power cord is attached to the receptacle on the right side of the device.
2. Plug the power cord into a Hospital Grade receptacle and switch on the power switch located in the power control module above the power cord. The light in the power switch will be green when the power is on.
3. The fan will start when the power switch is turned on.
4. Unlock the system door with the key supplied and push the key button/latch to open the door. Make sure the key is kept in a safe place.
5. On the upper right of the system is an indicator light to notify operator when it is necessary to change the HEPA filter.
6. When filter is within operating parameters (11-months) the indicator light will turn off.
7. Once eleven (11) months running time has elapsed, the indicator light will turn yellow indicating it is time to change the filter. See "*Filter Maintenance*" section for instructions on replacing the filter.


Probe Hanging: Vaginal, Rectal & General Purpose Probes Only

WARNING: Follow the probe manufacturer's instructions for storage and handling.

CAUTION: When installing the Vaginal/Rectal/General Purpose probes, ensure at all times that the connector housing and cable are kept away from the transducer to avoid cross-contamination.

1. Open the system door by depressing the key button/latch.
2. While holding the transducer in an upright position, guide it up through the circular slot in the upper shelf. Rotate the transducer distal tip/lens so that it is facing outward.
3. Guide the cable through the slots in the middle and lower shelves. Rest the transducer securely on the middle shelf.
4. Insert the connector housing into the circular slot on the lower shelf, with the cable looped and hanging in a downward position behind the connector.
5. Repeat steps 2 through 4, as needed.
6. Close door.

NOTE: Probe manufacturers will have the connector housing orientated in different positions. Ensure that the connector housing is secured safely on the lower shelf before moving to the next step.

FOR GENERAL PURPOSE PROBES:

1. Open the system door by depressing the key button/latch.
2. While holding the transducer, guide the cable through the slots in the middle and lower shelves. Rest the transducer securely on the middle shelf.
3. Insert the connector housing into the circular slot on the lower shelf, with the cable looped and hanging in a downward position behind the connector.

NOTE: Probe manufacturers will have the connector housing orientated in different positions. Ensure that the connector housing is secured safely on the lower shelf.

4. Repeat Steps 2 and 3, as needed.


Probe Hanging: TEE Probes Only

WARNING: Follow the probe manufacturer's instructions for storage and handling.

CAUTION: When installing the TEE probe, ensure at all times that the connector housing and cable are kept away from the shaft and transducer lens to avoid cross-contamination.

CAUTION: Ensure that the distal tip and transducer lens are stored securely away from the system door to avoid hitting the system door when closed.

1. Open the system door by depressing the key button/latch.
2. Hold the TEE Probe above the upper shelf while resting the shaft into the slot on the middle shelf, and resting the distal tip/transducer lens into the slot on the lower shelf.
3. Lower the control housing strain relief into one of the circular slots located on the upper shelf.
4. Take the cable and loop it over the strain relief located above the upper shelf.
5. Insert the cable into the slots located on the upper and middle shelf.
6. Guide the cable through the slot on the middle shelf, and place the connector housing on the upper shelf on the left side of the TEE probe.
7. Repeat Steps 2 through 6, as needed.
8. Close the door.

Cleaning and Disinfecting the System

WARNING: Users of this product have an obligation and responsibility to provide the highest degree of infection control to patients, co-workers and themselves. To avoid cross-contamination, follow infection control policies established by your facility.

CAUTION: Disconnect power when cleaning the system.

CAUTION: Do not use solvents or isopropyl alcohol to clean the acrylic door. Doing so may damage the doors. Use a mild detergent solution on a soft cleaning cloth.

CAUTION: Do not wipe surfaces with hard or abrasive materials. Doing so will scratch the surfaces.

1. To clean and disinfect system door and acrylic windows, use a soft cleaning cloth and mild detergent solution.

2. To clean and disinfect system, use a surface disinfectant containing up to 70% isopropyl alcohol, such as CiDehol® 70, authorized by your infection control department.
3. Make sure that the system is completely dry before use.
4. The system is able to be cleaned and disinfected weekly. Follow your facility's guidelines for frequency of cleaning and disinfecting the system.

Routine Maintenance

This system requires minimal routine maintenance and inspection.

Daily System Inspection

1. Visually inspect electrical connections prior to operation.
2. Check area around the device to ensure safe operation.
3. Check fan operation by placing hand over the outlet on the top of the inside of the cabinet and feeling for air movement.
4. Check fan operation by monitoring for unusual or harsh noises.
5. Confirm filter replacement indicator is off. If the indicator is yellow, the HEPA filter needs to be replaced. See "Filter Maintenance" section for instructions on replacing the filter.

Filter Maintenance

Filter life is monitored by a timer; replacement is signaled after the filter has been used for eleven (11) months. The timer is only active when the fan is on. If the Storage Cabinet were to lose power, the timer inside will reset to the last completed hour. Call your CIVCO Medical representative or local distributor to order replacement filters.

CAUTION: Filter life is dependent upon frequency and conditions of use. As a result, the filter may need to be replaced prior to the 11 month guideline.

CAUTION: Disconnect power when replacing the filter.

NOTE: The HEPA filter is fiberglass based and does not have a shelf life. The filter will perform as required if properly sealed and stored before use.

1. Turn the main power switch off and disconnect the power cord from the electrical outlet.
2. Due to the height of the device, a stool or ladder is recommended to access the filter.
3. There are two latches on the top panel of the device. To open the top and access the filter, turn the latches 90 degrees. Use the latches on the top panel to lift the top panel to a vertical position.
4. Remove the used filter by pulling it up.
5. Put the used filter into a plastic bag and dispose of it with other non-hazardous waste. **This does not need to be red tagged.**

-
6. Install the new HEPA filter with the gasket and arrow pointing down.

NOTE: HEPA filters are very fragile. When handling the HEPA filter, use the outer metal frame to grip the filter when changing or installing. Do not damage the filter material. Handle the filter on the outer metal frame.

7. Once the filter is in place, swing the top panel down, push down and secure by closing the left and right latches previously opened.
8. Pull up on the latches on the top of the device to make sure latches are secure.
9. Plug the device back into the wall outlet and turn on the power switch.

NOTE: The filter change indicator light will still be yellow at this point and will need to be reset.

10. To reset the filter change indicator light, located on the upper side of the system behind the power inlet module, hold the button down for 5 seconds until the maintenance indicator light briefly flashes then turns off. This will reset the filter change indicator light for another 11 months.

Technical Support

CIVCO Medical Solutions, 750 Lois Lane, Coralville, IA 52241 USA
Tel: 319.248.6757 / 800.445.6741
WWW.CIVCO.COM

Registration Card

Please complete and return the attached Warranty Card to CIVCO Medical Solutions. Call 1-800-445-6741 to purchase additional filters.

Appendix A: Troubleshooting

If the CIVCO Storage System fan is not running, check the following.

1. Ensure that the device is plugged into a live wall outlet.
2. Make sure the On/Off switch is on. The switch will be illuminated when on.
3. In the event that the system is still not working, you will need to check the fuses, which are externally mounted in the power entry module with the On/Off switch and power cord receptacle.
4. The fuses are located in a rectangular receptacle between the switch and the power cord. Remove the fuse receptacle by squeezing the two tabs on the fuse receptacle and pulling outward. If fuses are blown, replace fuse.
5. If the system is still not working, call CIVCO Medical Solutions at 800.445.6741 or call your local distributor.

NOTE: If the filter change indicator light is yellow, this indicates it is time to change the HEPA filter. Refer to "*Filter Maintenance*" section for instructions on replacing the filter.

Appendix B: Site Requirements

Actual Size

- Width: 24.0" 60.96 cm
- Depth: 10.0" 25.40 cm
- Height: 74.0" 187.96 cm
- Weight: 100 lbs 45.36 kg

Minimum Space Required

- Width: 30.0" 76.20 cm
- Depth: 36.0" 91.44 cm
- Height: 84.0" 213.36 cm

Electrical Requirements


- Service: 100-240 VAC 50/60Hz, 15A Hospital Grade GFCI receptacle
- Fuse: 1.0 Amp 250V 5mm x 20mm FSF Panel Mounted
- Current Draw: 0.5 Amp
- Power Input: 8.3 Watts

Environmental Conditions

- Indoor use only
- Altitude up to 3,000 meters
- Temperature range: 5°C to 40°C
- Maximum relative humidity: 80% for temperatures up to 31°C, decreasing linearly to 50% relative humidity at 40°C
- Main supply voltage fluctuations not to exceed +/- 10% of the nominal voltage
- Insulation over voltage category II
- Pollution degree 2

Mounting and Installation Diagrams

STORAGE CABINET SCHEMATIC

Appendix C: Product Specifications

Product Consumables

CIVCO Part #	Description
610-1328	Replacement Filter

Airflow Specifications

Volume of Filtered Air: 20 Air Exchanges per hour (ACH)

Fan Speed: 31 CFM (53 m³/h)

Noise Level Specifications

Sound Pressure Level: <45 dBA


Cabinet Materials

Part Description	Material
Interior of Unit	Powder Coated Sheet Metal
Exterior of Unit	Powder Coated Sheet Metal
Rivets	Zinc-Plated Steel
Door Handle (with rubber stopper)	Zinc-Plated Steel, Powder Coated Zinc, Nickel-Plated Brass, Neoprene
Hinges	Acetal, Polycarbonate
Filter Indicator Light Lens	Polycarbonate
Power Connector / Switch	Nylon
Lift Handles	ABS
Shelving	ABS
Window	Acrylic
Nuts and Screws	Stainless Steel
Door Bumpers	Polyurethane
Shelving Screw Caps	Plastic
Latches	Nickel-Plated Zinc


Shelving Specifications

Storage System for VR and General Purpose Probes


Top Shelf

Middle Shelf


Bottom Shelf


Storage System for TEE probes



Top Shelf

Middle Shelf

Bottom Shelf

Appendix D: EMC Compliance Information

The device has been tested and determined to be in compliance with the standards or test levels noted in Table 1 below with no deviations from the test levels noted.

No errors were produced as a result of EMC testing. However, consideration should be given to the EM environment the product is operated in. Portable RF communications equipment (including peripherals such as antenna cables and external antennas) should be used no closer than 30cm (12 inches) to any part of the system, including cables specified by the manufacturer. Degradation of the performance of this equipment could result.

Refer to AAMI TIR 28 for guidance in management of the EM environment and management of medical devices for EMC, including assessment of the EM environment, investigation and reporting of EMI problems.

Phenomenon	Test Levels
Radiated Emissions	CISPR 11 Class A
Conducted Emissions	CISPR 11 Class A
Harmonics	Class A
Voltage Fluctuations and Flicker	Compliant
Electrostatic Discharge (ESD)	±8 kV Contact ±2 kV, ±4 kV, ±8kV, ±15 kV air
Radiated Immunity (including Proximity Fields)	3 V/m 80 MHz to 2.7 GHz 80% AM at 1kHz
Immunity to RF Wireless Communication Equipment (frequencies)	IEC 60601-1-2 Table 9
Electrical Fast Transients/Bursts	±2 kV 100 kHz repetition frequency
Surges Line to Line	±0.5 kV, ±1 kV
Surges Line to Ground	±0.5 kV, ±1 kV, ±2 kV
Conducted Immunity	3 V 0.15 MHz - 80 MHz 6 V in ISM bands between 0.15 MHz and 80 MHz 80% AM at 1 kHz
Magnetic Immunity	30 A/m 50 Hz or 60 Hz
Proximity to Magnetic fields immunity	8 A/m @ 30 kHz (CW) 65 A/m @ 134.2 kHz (50% Pulse at 2.1 kHz) 7.5 A/m @ 13.56 MHz (50% Pulse at 50 kHz)
Voltage Dip	0% Ut: .5 cycle At 0°, 45°, 90°, 135°, 180°, 225°, 270°, 315°
	0% Ut: 1 cycle 70% Ut: 25/30 cycles Single phase at 0°
Voltage Interruptions	0% Ut: 250/300 cycles

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