


Imaging Table with Fowler Back and Drop Window



MODEL 80072

- Two-piece top with pneumatic adjustable backrest and window drop
- Trendelenburg/Reverse Trendelenburg plus or minus 15°
- Easy-to-use hand control
- All-welded base and top frame
- Unobstructed base design
- Dual column lift system adds rigidity
- Paper dispenser included
- Industry leading 5-year, limited warranty
- 34" wide model available
- Model 80072: 500 lbs load capacity under normal use
- Model 80072-X: 600 lbs load capacity under normal use

California residents please be advised of the following, pursuant to Proposition 65:



WARNING: This product can expose you to chemicals, including chromium, which is known to the State of California to cause cancer or reproductive toxicity. For more information go to www.P65Warnings.ca.gov.



Trendelenburg



Reverse Trendelenburg



Optional Side Rails & Casters

80072 SPECIFICATIONS

	length	width	height range
overall	72"	27"	21 ¹ / ₂ " – 37 ¹ / ₂ "
backrest	38"	27"	
window	10"	10"	

80072-X Extra-Wide 34" width model



Upholstery Colors



OPTIONS

Safety Straps Set

- 2" x 92" strap is secured with quick hook and loop Velcro® fasteners
- 2 safety straps per set

080	length	width
	92"	2"
	233.67 cm	5 cm

Safety Rails (pair)

- 3 height positions
- Plunger height adjustment

098

Power Backrest

- Smooth, power control raises to 80° and lowers to flat

011

Locking Casters

- Heavy duty, rubber wheel, easy-roll casters

087

Paper Cutter

- Clear plastic with elastic straps and buckle

040

