

# SONOGRAPHER'S CHAIR

Sonography chairs provide adjustments for seat height, back/arm height, and backward tilt tension, along with an integrated seat slider. These adjustments, along with the ergonomically designed support surfaces, allow sonographers to find the least stressful position for a wide variety of scans.



**Capisco 8106 Model**  
\*Available with optional foot ring



**Capisco 8010 Model**  
\*Shown with optional foot ring.

## Features

- Made for greater freedom of movement, variation and natural sitting positions
- Gives you an upright sitting posture and a natural curvature in your lower back
- Designed to be used facing forward, backward or side saddle
- Seat height, seat depth, back height, backward tilt tension and lockable tilt

Specifications	8106	8010
Height Medium	19" - 26"	18" - 24.5"
Height Tall	22.5" - 32"	21" - 30"
Chair Width	23.5"	22"
Chair Weight	31 lbs.	24.5 lbs.
Weight Rating	300 lbs.	300 lbs.
Foot Ring	Optional	Optional

717.759.3146 [www.oakworksmed.com](http://www.oakworksmed.com) email: [medsales@oakworks.com](mailto:medsales@oakworks.com)