

# PRODUCT DATASHEET

Model: CS2

Description: Clear 2 Step Imaging Footstool

Revision: October 2020



## Product Attributes

- Ideal for lateral and AP weight bearing exams
- 3 separate slot positions for image receptor
- Safety handles for patient support
- Durable, non-slip surfaces
- 16 in. high platform
- Can accommodate a 17 x 17 in. cassette size imaging receptor or smaller
- Tilt-back casters for easy mobility
- Rated up to 500lb.

## Material Composition

Material	%
Aluminum	58%
Polycarbonate	15%
HDPE	17%
Steel	8%
Vinyl	1%

## DIMENSIONS

Base Outside Dimensions: 26.25 x 32.5in. (66.68 x 82.55cm)

Handrail Height (from floor): 52.13in. (132.4cm)

Platform Height: 16.5in (41.9cm)

Platform Dimensions: 24.88 x 21in. (63.2 x 53.34cm)

Lateral Slot Dimensions: 19 x 1in. (48.26 x 2.54cm)

AP Tray Dimensions: 22.88 x 1.13in. (58.12 x 2.87cm)

Overall Weight: 62lb. (28kg)

## STATEMENT OF CONFORMANCE

Reina Imaging hereby states that the data, specifications, and descriptions contained herein this document are accurate and representative of the models being described.

## REGISTRATIONS AND COMPLIANCE

All models listed above are made in the USA and registered to:

Reina Imaging

6107 Lou Ave.

Crystal Lake, IL 60014 USA