

MOBILE LEADED BARRIER WITH 12" X 10" WINDOW

PRODUCT INFORMATION

COPYRIGHT ©2025. ALL RIGHTS RESERVED



PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW




The Mobile Leaded Barrier with 12" x 10" Window is a portable, full-body radiation barrier designed to protect against secondary radiation. It features 2.0mm lead equivalent opaque panels of solid, leaded x-ray protection. Its 12" W x 10" H lead equivalency window is a viewing panel with 2.0mm lead equivalent Lead glass. Available in eight different sizes (24" x 72", 24" x 84", 30" x 72", 30" x 84", 36" x 72", 36" x 84", 45" x 72", and 45" x 84"), this Phillips Safety leaded barrier travels easily on its four-caster locking base steel legs. The window centers are 60" from the bottom of the panel and it is recommended that two people move the barrier together at all times. Please note that 84" high barriers may not fit through all doorways, therefore make sure you have selected the correct size.

SPECIFICATION

MATERIAL: Lead Glass with opaque panel solid lead	PROTECTION: 2.0mm Pb LE opaque panels
FEATURES: Open base design, portable, multiple sizes available, heavy duty, 4-castor base	SIZES AVAILABLE: 24" x 72" 24" x 84" 30" x 72" 30" x 84" 36" x 72" 36" x 84" 45" x 72" 45" x 84"
	WINDOW SIZE: 12" x 10"

PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW

Mobile Leaded Barrier
with 12" x 10" WINDOW





LB-2472

SCALE: 1/2" = 1'

PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW

**Mobile Leaded Barrier
with 12" x 10" WINDOW**




LB-2484

SCALE: 1/2" = 1'

PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW


Mobile Leaded Barrier
with 12" x 10" WINDOW




LB-3072

SCALE: 1/2" = 1'

PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW






PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW



PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW

**Mobile Leaded Barrier
with 12" x 10" WINDOW**




LB-3684

SCALE: 1/2" = 1'


PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW

Mobile Leaded Barrier
with 12" x 10" WINDOW



Top



Front



Side



LB-4572

SCALE: 1/2" = 1'

PRODUCT INFORMATION
MOBILE LEADED BARRIER WITH 12" X 10" WINDOW

